

R14374, R14689

FEATURES

- ●High gain 1.0 × 10⁷
- ϕ 80 mm Hemisphere...R14374

APPLICATIONS

High energy physics



Left: R14374, Right: R14689

SPECIFICATIONS

GENERAL

GENTENAL						
	Parameter	R14374	R14689	Unit		
Spectral response		300 to	nm			
Wavelength of maximum response		42	nm			
Window material		Borosilica				
Photocathode	Material	Biall				
	Minimum effective area	φ72	φ81	mm		
Dynode	Structure	Circular and li	_			
	Number of stages	10				
Base		JEDEC No	_			
Operating ambient temperature		-30 to	°C			
Storage temperature		-30 to	°C			
Suitable socket		F678-14W (So	_			

MAXIMUM RATINGS (Absolute maximum values)

	Parameter	R14374	R14689	Unit
Supply voltage	Between anode and cathode	1500		
	Between anode and last dynode	300		
Average anode current		0.1		

CHARACTERISTICS (Typ.) (at 25 °C)

	Parameter	R14374	R14689	Unit		
Cathode sensitivity	Luminous (2856 K)	90				
	Radiant at 420 nm	90				
	Blue sensitivity index (CS 5-58)	11.0				
	Quantum efficiency at 380 nm	27.5				
Anode sensitivity	Luminous (2856 K)	900				
	Radiant at 420 nm	9.0 >	A/W			
Gain		1.0 >				
Anode dark current (A	After 30 minute storage in darkness)	5	nA			
Time response	Anode pulse rise time	2.9	2.9	ns		
	Electron transit time	35	36	ns		
	Transit time spread (FWHM)	1.3	1.5	ns		

NOTE: Anode characteristics are measured with a voltage distribution ratio and supply voltage shown below.

VOLTAGE DISTRIBUTION RATIO AND SUPPLY VOLTAGE

Electrodes	K	Dy1	Dy2	Dy3	Dy4	Dy5	Dy6	Dy7	Dy8	Dy9	Dy10	Р
Ratio	(3	1	1	1	1	1	1	1	1	1	1

Supply voltage: 1250 V, K: Cathode, Dy: Dynode, P: Anode

Subject to local technical requirements and regulations, availability of products included in this promotional material may vary. Please consult with our sales office. Information furnished by HAMAMATSU is believed to be reliable. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications are subject to change without notice. No patent rights are granted to any of the circuits described herein. ©2020 Hamamatsu Photonics K.K.

PHOTOMULTIPLIER TUBE R14374, R14689

Figure 1: Typical spectral response




Figure 2: Typical gain




Figure 3: Transit time spread (FWHM)







Figure 4: Dimensional outline (Unit: mm)

●R14374



●R14689





TPMHA0635EC

HAMAMATSU PHOTONICS K.K. www.hamamatsu.com

Electron Tube Division

314-5, Shimokanzo, Iwata City, Shizuoka Pref., 438-0193, Japan, Telephone: (81)539/62-5248, Fax: (81)539/62-2205

U.S.A.: Hamamatsu Corporation: 360 Foothill Road, Bridgewater, NJ 08807, U.S.A., Telephone: (1908-231-0960, Fax: (1908-231-1218 E-mail: usa@hamamatsu.com

Germany: Hamamatsu Photonics Deutschland GmbH:. Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-265-8 E-mail: info@hamamatsu.de

France: Hamamatsu Photonics France S.A.R.L.: 19, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 10, Fax: (33)1 69 53 71 10 E-mail: info@hamamatsu.fr

United Kingdom: Hamamatsu Photonics UK Limited'z Howard Court,10 Tewin Road, Wellyn Road, Wellyn Tor 294888, Fax: (44)1707-294888, Fax: (44)1707-324888, Fax: (44)1707-